



The University of Texas at Austin
**Aerospace Engineering
and Engineering Mechanics**
Cockrell School of Engineering

ASE 375 Electromechanical Systems
Section 14115

Monday: 3:00 - 6:00 pm

Report 7: Shake Testing

Andrew Doty, Andres Suniaga, Dennis Hom
Due Date: 04/01/2024

Contents

| | | |
|----------|-----------------------------------|-----------|
| 1 | Introduction | 2 |
| 2 | Equipment | 2 |
| 3 | Procedure | 4 |
| 3.1 | Load Cell Calibration | 4 |
| 3.2 | Shake Testing | 5 |
| 4 | Data Processing | 6 |
| 4.1 | Variables and Equations | 6 |
| 5 | Results and Analysis | 7 |
| 5.1 | Part 1 | 7 |
| 5.2 | Part 2 | 8 |
| 6 | Conclusion | 9 |
| A | Datasheets | 11 |

1 Introduction

In this experiment we analyze the frequency response of a scaled-down wing model using a load cell, MEMS accelerometer, and electromagnetic shaker. The goal of this lab is to learn how a load cell, electromagnetic shaker, and instrumentation amplifier operate by performing a shake test on the scaled-down wing.

The shaker will vibrate the scaled wing model at certain frequencies and the load cell will be measuring the force while the accelerometer will be measuring the motion of the oscillation. From this we will obtain a frequency response function for the wing. This will provide useful information of the wing's structural properties.

2 Equipment

Measurement devices and hardware used in this lab include:

- Built-Up Wing Model:

Scaled-down wing model used for shake testing.




Figure 1: Built-Up Wing Model Setup, for Part 2 of this Experiment

- Strain Gage Load Cell [5]:

Sensor that measures force as an electrical output. Within the load cell there are strain gages set up in a Wheatstone bridge configuration making the load cell's operating principle the variable resistance caused by application of force or deformation which is then measured as an electric signal. The load cell is calibrated in the first part of the lab as shown in Figure 2 and then placed on a rod atop the shaker to attach to the wing as shown in Figure 1.

- 'Rare Earth' Magnets:

Small Magnets used for secure placement of load cell under the wing platform.

- Brass Slotted Weights with hanger and L-bracket:

L-bracket used for placement of load cell on beam. Weights hanged on L-bracket to calibrate the load cell as shown below:




Figure 2: Placement for Calibration of the Load Cell, for Part 1 of this Experiment

- MEMS Accelerometer (ADXL335) [2]

The ADXL335 is a small 3-Axis MEMS Accelerometer which operates based on differential capacitance. An applied force or acceleration will change the capacitance within the MEMS accelerometer. In this experiment, the MEMS accelerometer is placed as shown in Figure 1, using only the Y-axis for our experiment. Operational Voltage = 3.3V

- Mechanical Wave Driver [4]:

Electromagnetic shaker with driver arm capable of oscillating/vibrating at variable amplitudes and any frequency between 0.1 Hz to 5 kHz. In this experiment it is used to shake the wing model as shown in Figure 1.

- DAQ, NI-9215 Voltage Input Module [1], NI-9263 Analog Output Module [6], and LabVIEW:

Data Acquisition System used to process sample measurements into digital data for the computer to read.

NI-9215 is an analog input module used to measure the output voltage signals of sensors and send it through the DAQ system. The NI-9215 module will measure the output of the instrumentation amplifier [3].

The NI-9263 output module will connect to the mechanical wave driver to supply it with sinusoidal waves.

LabVIEW used to model these output voltages read from the DAQ of the accelerometer and load cell measurements.

- Solderless Breadboard, Jumper Wires, and AD623 Instrumentation Amplifier [3] :

Used to make connections to the input analog modules and to construct circuits. In this lab we connect the accelerometer inputs to the breadboard for power, ground, and signal for the measurement axis.

We also make a circuit for the load cell as shown below:




Figure 3: Load Cell Amplifier Circuit

3 Procedure

3.1 Load Cell Calibration

Before beginning the calibration, ensure the block diagram to gather the data on LabVIEW is setup as shown in the example below:




Figure 4: LabVIEW Block Diagram Setup Example

- (DAQ) Data → Signals (Collector),
- (Collector) Collected Signals → Signals (Statistics),
- (Collector) Collected Signals → Waveform Graph,
- Arithmetic Mean → Mean,
- Standard Deviation → STD.

This LabVIEW setup will plot the load cell compression and expansion as an electrical output over time as well as calculate the average voltage, V_{OUT} , and standard deviation. These are the values we need to mark down for calibration.

NOTE: Disregard "Mems(Y)" name over the waveform plot. Rename as "Load Cell".

Calibrating the load cell:

1. First, with the appropriate screws and L-bracket, screw in the top of the load cell to a beam with the L-bracket screwed just underneath the load cell as shown in Figure 2. Weights will be hung on the L-bracket to calibrate the load cell.
2. Measure the gain setting resistor, R_G as shown in Figure 3.
NOTE: This step can be done before or after calibrating of the load cell.
3. Run the LabVIEW model. Load sequence of weights (in grams) to hang on the L-bracket,

$$\mathbf{w} = [0, 50, 70, 90, 110, 130, 150, 170, 190, 210, 230, 250, 270]$$

Zero means no weight on the bracket and 50 grams is the weight of the hangar by itself.

4. Mark down measurements of the average electrical output (mean) and the standard deviation.
5. This data can be found in Data Processing

3.2 Shake Testing

Block Diagram: Open "sineswept.vi" for LabVIEW. This will allow the software to send sinusoidal waves through a range of frequencies from 5 Hz to 100 Hz to the shaker, which we will use for this experiment. Use 1 kHz sampling frequency.

Shake Test:

1. Connect the MEMS accelerometer and load cell to the circuit board (amp circuit for load cell shown in Figure 3) and respective NI-9215 input modules. Ensure Shaker is connected to the NI-9263 output module to receive sinusoidal waves of different frequencies. Setup the wing model for shake testing as shown in Figure 1.
2. Issues to prevent before beginning shake test:
 - Load Cell Slippage off the wing model
 - Bent Rod joining the load cell and mechanical wave driver
 - Slanted or Loose Accelerometer
3. Adjust the amplitude for the wave driver as needed. For this experiment, our amplitude was placed on the lower end of the knob values.
4. Perform shake test by running the LabVIEW Sine Swept model. The wing will start shaking and frequency will increase until 100 Hz is reached.
5. Restart shake test if load cell slips off or accelerometer comes loose.
6. Once the shake test concludes, save the data which will be written to a file.
7. Next up is Data Processing to find the Frequency Response function of the wing.

4 Data Processing

4.1 Variables and Equations

Gain of the AD623 Instrumentation Amplifier:

$$G = 1 + \frac{100 \text{ k}\Omega}{R_G} \quad (1)$$

Variables:

- R_G : Resistance of Gain setting resistor

Differential Voltage, Output Voltage (read by DAQ) is amplified by Gain G :

$$V_{\text{OUT}} = GV_0 \quad (2)$$

Variables:

- V_{OUT} : Output voltage read by DAQ
- V_0 : Input voltage going into Instrumentation Amp

Harmonic Motion Equations:

$$x(t) = A \sin(\omega t) \quad (3)$$

$$\ddot{x}(t) = -A\omega^2 \sin(\omega t) \quad (4)$$

With these harmonic motion equations we can make the estimation:

$$|x| = A, \quad |\ddot{x}| = A\omega^2, \quad \implies |x| = \frac{|\ddot{x}|}{\omega^2}$$

Variables:

| Mass | Weight (N) | V (mean) | V (STD) | V (norm) |
|------|------------|----------|---------|----------|
| 0 | 0 | -0.00172 | 0.0005 | 0 |
| 50 | 0.491 | 0.00282 | 0.0006 | 0.00454 |
| 90 | 0.883 | 0.00647 | 0.0005 | 0.00819 |
| 130 | 1.275 | 0.0102 | 0.0005 | 0.0119 |
| 170 | 1.667 | 0.0141 | 0.0005 | 0.0158 |
| 210 | 2.06 | 0.0176 | 0.0004 | 0.0193 |
| 250 | 2.453 | 0.0215 | 0.0005 | 0.0232 |

Table 1: Load Cell Calibration

- $x(t)$: Displacement as a function of time
- $\ddot{x}(t)$: Acceleration as a function of time
- A : Amplitude
- ω : Frequency

5 Results and Analysis

5.1 Part 1




Figure 5: Load Cell Calibration Graph

This leads to a Calibration Factor of

$$C = 0.0095 \frac{V}{N} \quad (5)$$

Through the use of a multimeter, R_G is measured as 2.699k Ω . To calculate the gain of the amplifier, **Eq. 1** and **Eq. 2** are used. G is calculated to be 38.051.

5.2 Part 2




Figure 6: Tip Displacement vs. Frequency




Figure 7: Force vs. Frequency




Figure 8: Acceleration vs. Frequency

For Figures 6,7, and 8, the behavior changes between 50 and 65Hz. For the Acceleration graph, the data is near-zero except for within this range of frequencies, where there is a sharp jump. Since the tip displacement is calculated from the acceleration, it follows the same trend. For Figure 7, the force is near constant until it starts to dip within this range of frequencies, reaching a minimum at 55Hz.




Figure 9: Frequency Response Function

The Frequency Response function is the tip bending displacement vs. input force. It is taken by dividing the magnitude of the tip displacement by the input force and then plotting it for each frequency tested. This is shown in Figure 9, where the peak of this graph is the natural frequency, which is 55 Hz. This is the same frequency where the force was at a minimum, and within the range where the acceleration and tip displacement were at their largest values.

6 Conclusion

In this lab experiment we learned how an electromagnetic shaker and load cell work in order to explore the dynamic response of the scaled-down wing. An AD623 instrumentation amplifier was used to amplify the output of the load cell. Thorough shake testing, a frequency response function was obtained for the wing. This experiment gave a natural frequency of 55 Hz. The previous lab had a calculated natural frequency of 59 Hz.

Appendix

A Datasheets

- [1] **NI-9215 Datasheet:**
<https://www.amc-systeme.de/files/pdf/ni-9215-amc.pdf>
- [2] **ADXL335 Accelerometer Datasheet:**
<https://www.analog.com/media/en/technical-documentation/data-sheets/adxl335.pdf>
- [3] **AD623 Instrumentation Amplifier Datasheet:**
<https://www.analog.com/media/en/technical-documentation/data-sheets/ad623.pdf>
- [4] **Mechanical Wave Driver (SF-9324) Specs:**
https://cdn.pasco.com/product_document/Mechanical-Wave-Driver-Manual-SF-9324.pdf
- [5] **Miniature Low-Profile Tension Link Load Cell Specs:**
<https://br.omega.com/omegaFiles/pressure/pdf/LCM703.pdf>
- [6] **NI-9263 Datasheet:**
<https://www.ni.com/docs/en-US/bundle/ni-9263-specs/page/specs.html>